
# Tyvek® IsoClean®

## Cat. III PROTECTION LEVEL


- Instructions for Use
- Gebrauchsanweisung
- Consignes d'utilisation
- Istruzioni per l'uso
- Instrucciones de uso
- Instruções de utilização
- Gebruiksaanwijsties
- Bruksanvisning
- Brugsanvisning

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ENGLISH

### INSTRUCTIONS FOR USE

**INSIDE LABEL MARKINGS** ① Trademark. ② Coverall manufacturer. ③ Tyvek® IsoClean® model identification – see table. ④ Processing and packaging identification – see table. ⑤ CE marking – Overall compliance with requirements for category III personal protective equipment according to European legislation. ISO 16602:2007/Amd.1:2012, type-test and quality assurance certificates were issued by SGS United Kingdom Ltd., Weston-super-Mare, BS22 6WA, UK, identified by the EC Notified Body number 0120. ⑥ Indicates compliance with European and ISO standards for chemical protective clothing. ⑦ Protection against particulate radioactive contamination according to EN 1073-2:2002. ⑧ EN 1073-2 clause 4.2 requires puncture resistance of class 2. This coverall meets class 1 only. ⑧ Full-body protection 'Types' achieved by this coverall defined by the European and ISO standards for Chemical Protective Clothing. ⑨ Wearer should read these instructions for use. ⑩ Sizing pictogram indicates body measurements. Check your body measurements and select the correct size. ⑪ Do not re-use. ⑫ Flammable material, keep away from fire. ⑬ This pictogram indicates protection against biological hazard.

THE FIVE CARE PICTOGRAMS INDICATE:

Do not wash.	Do not iron.	Do not machine dry.	Do not dry clean.	Do not bleach.

MODEL IDENTIFICATION DETAILS:

<b>IC183B</b>	IC183B is the model name for a sterile non-hooded protective coverall with bound seams, cuff, ankle and waist elastic.
---------------	--

PROCESS AND PACKAGING IDENTIFICATION DETAILS:

<b>DS</b>	clean processed, sterilized and double bagged.
-----------	--

PERFORMANCE OF STERILIZED TYVEK® AND TYVEK® ISOCLEAN® COVERALLS:

FABRIC PHYSICAL PROPERTIES	TEST METHOD	RESULT	UNIT	EN CLASS*
Abrasion resistance	EN 530 method 2 ISO 12947-2	> 10	Cycles	1/6
Flex cracking resistance	EN ISO 7854 method B	> 100000	Cycles	6/6
Trapezoidal Tear resistance	EN ISO 9073-4	> 10	N	1/6
Tensile strength	EN ISO 13934-1	> 30	N	1/6
Puncture resistance	EN 863	> 5	N	1/6
Burst resistance	EN ISO 13938-1	> 80	kPa	2/6

\*According to EN 14325:2004 and ISO 16602:2007/Amd.1:2012

RESISTANCE TO PENETRATION BY LIQUIDS (EN ISO 6530)

Chemical	Penetration index - EN Class*	Repellency index - EN Class*
Sulphuric acid (30%)	3/3	3/3
Sodium hydroxide (10%)	2/3	2/3

\*According to EN 14325:2004 and ISO 16602:2007/Amd.1:2012

FABRIC RESISTANCE TO PENETRATION OF INFECTIVE AGENTS

Test	Test method	EN Class*
Resistance to penetration by blood and body fluids using synthetic blood	ISO 16603 Procedure C	3/6
Resistance to penetration by blood-borne pathogens using Phi-X174 bacteriophage	ISO 16604 Procedure C	no classification
Resistance to penetration by contaminated liquids	EN ISO 22610	1/6
Resistance to penetration by biologically contaminated aerosols	ISO/DIS 22611	1/3
Resistance to penetration by contaminated solid particles	ISO 22612	1/3

\* According to EN 14126:2003

WHOLE SUIT TEST PERFORMANCE

Method	Test result	EN Class
Type 5: Particulate aerosol inward leakage test (EN 13982-2)	Pass*** $L_{inj} 82/90 \leq 30\%^{**}$ $L_8/10 \leq 15\%^{**}$	N/A
Protection factor according to EN 1073-2:2002	> 50	2/3***
Type 6: Low level spray test (EN ISO 17491-4:2008, method A)	Pass****	N/A
Seam strength (EN ISO 13935-2)	> 30	1/6*

N/A=Not Applicable \*According to EN 14325:2004 and ISO 16602:2007/Amd.1:2012 \*\* 82/90 means 91,1%  $L_{inj}$  values  $\leq 30\%$  and 8/10 means 80%  $L$  values  $\leq 15\%$

\*\*\*Test performed with taped cuffs, hood, ankles and zipper flap. \*\*\*\*Test performed with separate Tyvek® hood.

For further information, please contact your supplier or DuPont.

**TYPICAL AREAS OF USE:** DuPont™ Tyvek® IsoClean® coveralls are designed to help protect sensitive products and processes from contamination by people and to protect workers from some hazardous substances. They are typically used, depending on the chemical toxicity and exposure conditions, for protection against particles (Type 5), limited liquid splashes or low level sprays (Type 6).

**LIMITATIONS OF USE:** Do not re-use. Do not use the product if it has surpassed the expiry date. For sterilized coveralls if the packaging has been damaged and is no longer air-tight the product is no longer sterile. Do not re-sterilize the product. Tyvek® IsoClean® fabric and coveralls are not flame resistant and should not be used around heat, open flame, sparks or in potentially flammable or explosive environments. Tyvek® melts at about 135°C. These models do not comply with standard EN 1149-5 (surface resistance) and are unsuitable for use in explosive zones. Exposure to certain very fine particles, intensive liquid sprays and splashes of hazardous substances may require coveralls of higher mechanical strength and barrier protection than those offered by Tyvek® IsoClean® fabric and coveralls. Bound seams of Tyvek® IsoClean® do not offer a barrier to infective agents. For increased protection the user should select a coverall comprising seams that offer equivalent protection as the fabric (such as stitched and overlapped seams). For enhanced protection in certain applications, taping of cuffs, ankles, hood and zipper flap should be considered. Models with thumb holders should only be used with a double glove system, where the wearer puts the thumb holder over the under glove and the second glove should be worn over the coverall. Please ensure that you have chosen the Tyvek® IsoClean® model suitable for your job. The user shall be the sole judge for the correct combination of full body protective coveralls and ancillary equipment (gloves, footwear, respiratory protective equipment, etc.) and for how long a Tyvek® IsoClean® coverall can be worn on a specific job with respect to its protective performance, wear comfort or heat stress. DuPont shall not accept any responsibility whatsoever for improper use of Tyvek® IsoClean® coveralls.

**TYVEK® ISOCLEAN® MODEL SPECIFIC INFORMATION:**

MODEL	
IC183B DS	When using the unhooded Tyvek® IsoClean® coverall a separate Tyvek® IsoClean® hood, fully taped to the coverall should be worn to achieve the claimed whole suit protection levels. Ensure that the hood has a good fitting of facial opening, adequate shoulder coverage that should be worn under the coverall.

**PREPARING FOR USE:** It is the responsibility of the user to be trained in the donning, doffing, proper use, handling, storage, maintenance and disposal of Tyvek® IsoClean® coveralls. In the unlikely event of defects, do not wear the coverall.

**STORAGE:** The coveralls may be stored between 15 and 25 °C in the dark (cardboard box) with no UV light exposure. DuPont suggests that Tyvek® IsoClean® coveralls be used within 5 years, provided they are properly stored and pass a full visual inspection. High temperature, oxidizing gases, wet, cold, ultraviolet and ionizing radiation can significantly impact the long-term life of coveralls made by Tyvek® fabric. See expiry date on bag label.

**DISPOSAL:** This coverall can be recycled, incinerated or buried in a controlled landfill. Disposal restrictions depend upon the contamination incurred during use and are subject to national or local legislation.

For further information about the coverall and its barrier performance, please contact your Tyvek® IsoClean® supplier or visit [www.jp.dupont.com](http://www.jp.dupont.com).

The content of this instruction sheet was last verified by the notified body SGS in September 2015.

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**PROCESS AND PACKAGING IDENTIFICATION DETAILS:**

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Copyright © 2015 DuPont. All rights reserved. The DuPont Oval Logo, DuPont™, The miracles of science™ and all products denoted with ® or ™ are registered trademarks or trademarks of E. I. du Pont de Nemours and Company or its affiliates.